

ANALIZA BIAŁEK



100012944

## Laboratory of Analysis

**Client:** RIVER PEPTIDES  
**Order number:** 100012944  
**Name:** KLOW KL800117  
**Number of samples:** 1 vial

**Batch number:** KL800117  
**Receive date:** 2026-02-27  
**Start date:** 2026-03-03  
**Enda date:** 2026-03-09

### Sample identification

Substance	Method of taking a sample	Form
-----	Sample provided by customer	Blue powder in glass vial


### Results

Substance	Sample name	Content	Purity
BPC-157	KLOW KL800117	11,07mg	99%
TB-500	KLOW KL800117	10,87mg	99%
GHK-Cu	KLOW KL800117	51,18mg	99%
KPV	KLOW KL800117	10,81mg	99%

### Research comment

The test sample was applied in the amount of 3  $\mu$ l to the BioBasic 8 Thermo Scientific column. During the analysis with the method a linear gradient in the solutions of H<sub>2</sub>O and Acetonitrile with the addition of 0.1% trifluoroacetic acid, a peaks of the sought substances appeared. Its concentration, when compared to the reference chromatogram, show the table above. The measurement error is + - 0,05mg. The actual state of the sample sent before the test is presented on the photo below

### Photo of the test sample




### Password

K7B4FVPQ



**Analiza Białek sp. z o.o.**  
**KRS:** 0000966168  
**NIP:** 8992922970  
**REGON:** 521751875

**Headquarters and Laboratory:**  
ul. Duńska 9, 54-427 Wrocław  
Budynek DELTA

**Lab technician:**  
Analysis and research specialist  
mgr Arkadiusz Zajęc




## Chromatogram



**Password**

K7B4FVPQ



## Method

Time	H2O + 0,1% TFA	ACN + 0,1% TFA
0min	95%	5%
2min	95%	5%
12min	5%	95%
15min	5%	95%

### Password

K7B4FVPQ