

ANALIZA BIAŁEK



100017339

Laboratory of Analysis

Client: River Peptides
Order number: 100017339
Name: KLOW
Number of samples: 1 vial

Batch number: KL800717
Receive date: 2026-08-17
Start date: 2026-08-21
End date: 2026-08-27

Sample identification

Substance	Method of taking a sample	Form
-----	Sample provided by customer	Blue powder in glass vial


Results

Substance	Sample name	Content	Purity
BPC-157	KLOW 80MG KL800717	10,82mg	99%
TB-500	KLOW 80MG KL800717	10,72mg	99%
GHK-Cu	KLOW 80MG KL800717	58,93mg	99%
KPV	KLOW 80MG KL800717	11,24mg	99%

Research comment

The test sample was applied in the amount of 3 μ l to the BioBasic 8 Thermo Scientific column. During the course of the analysis with the method a linear gradient in the solutions of H₂O and Acetonitrile with the addition of 0.1% trifluoroacetic acid, a peaks of the sought substances appeared. Its concentration, when compared to the reference chromatogram shows the table above. Height the measurement error is + - 0,05 mg. The actual state of the sample sent before the test is presented on the photo below

Photo of the test sample




Password

S17U3K9G

Analiza Białek sp. z o.o.
KRS: 0000966168
NIP: 8992922970
REGON: 521751875

Headquarters and Laboratory:
ul. Duńska 9, 54-427 Wrocław
Budynek DELTA

Lab technician:
Analysis and research specialist
mgr Arkadiusz Zajac





ANALIZA BIAŁEK



100017339

Chromatogram



Password

S17U3K9G

Analiza Białek sp. z o.o.

KRS: 0000966168

NIP: 8992922970


REGON: 521751875

Headquarters and Laboratory:

ul. Duńska 9, 54-427 Wrocław
Budynek DELTA

Lab technician:

Analysis and research specialist
mgr Arkadiusz Zająć



ANALIZA BIAŁEK



100017339

Method

Time	H ₂ O + 0,1% TFA	ACN + 0,1% TFA
0min	95%	5%
2min	95%	5%
12min	5%	95%
15min	5%	95%

Password

S17U3K9G

Analiza Białek sp. z o.o.

KRS: 0000966168

NIP: 8992922970

REGON: 521751875

Headquarters and Laboratory:

ul. Duńska 9, 54-427 Wrocław
Budynek DELTA

Lab technician:

Analysis and research specialist
mgr Arkadiusz Zajac